

COMPATIBILITY OF

AHN myTip LFT Low Retention Filter Tips

AHN myTip LFT Low Retention Filter Tips Compatibility with AHN Single Channel Pipettes

Filter Tips Compatibility
AHN Single Channel Pipettes

Pipettor		AHN																													
Brand		pipet4u® pro										pipet4u® Revolution																			
Model		variable/single channel						fix/single channel					variable/ electronic single channel																		
Volume Range / Cat. No.		0.1-2.5 µl / 8-100-00-9	0.5-10 µl / 8-101-00-9	2-20 µl / 8-102-00-9	5-50 µl / 8-103-00-9	10-100 µl / 8-104-00-9	20-200 µl / 8-105-00-9	50-200 µl / 8-105-01-9	100-1000 µl / 8-106-00-9	200-1000 µl / 8-106-01-9	500-5000 µl / 8-107-00-9	2000-10000 µl / 8-108-00-9	5 µl / 8-001-00-9	10 µl / 8-002-00-9	20 µl / 8-003-00-9	25 µl / 8-004-00-9	50 µl / 8-007-00-9	100 µl / 8-008-00-9	200 µl / 8-009-00-9	250 µl / 8-010-00-9	500 µl / 8-011-00-9	1000 µl / 8-012-00-9	2000 µl / 8-013-00-9	2500 µl / 8-014-00-9	5000 µl / 8-015-00-9	10000 µl / 8-016-00-9	0.2-10 µl / 8-500-00-0	5-100 µl / 8-501-00-0	10-200 µl / 8-502-00-0	50-1000 µl / 8-503-00-0	100-5000 µl / 8-505-00-0
Volume		Premium-Line Filter Tips																													
Cat.No.																															
0.1-10 µl	5-001-C5-0		■	■																						■					
0.1-10 µl XL	5-061-C5-0	■	■	■																						■					
0.5-10/20 µl	5-051-C5-0		■																												
1-20/30 µl	5-116-C5-0				■	■	■	■				■	■	■	■												■	■			
1-50 µl	5-111-C5-0				■	■	■	■				■	■	■	■	■											■	■			
1-100 µl	5-112-C5-0				■	■	■	■				■	■	■	■	■											■	■			
1-200 µl	5-113-C5-0					■	■	■								■	■										■	■			
100-1000 µl long	5-201-C5-0									■	■						■	■	■											■	

Information correct at time of printing, Date: 29 May 17



metal-free

autoclavable
121° C | 15 min



AHN myTip LFT Low Retention Filter Tips Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN											
Brand		pipet4u [®] pro		pipet4u [®] Revolution									
Model		variable/multi channel		variable/multi channel									
		8 Channel	12 Channel	8 Channel	12 Channel								
Volume Range / Cat. No.		0.5-10 µl / 8-201-00-9	5-50 µl / 8-202-00-9	30-300 µl / 8-203-00-9	0.5-10 µl / 8-401-00-9	5-50 µl / 8-402-00-9	30-300 µl / 8-403-00-9	0.2-10 µl / 8-520-00-0	2.5-50 µl / 8-521-00-0	15-300 µl / 8-522-00-0	0.2-10 µl / 8-530-00-0	2.5-50 µl / 8-531-00-0	15-300 µl / 8-532-00-0
	Premium-Line												
Volume	Cat.No.												
0.1-10 µl	5-001-C5-0	■			■			■			■		
0.1-10 µl XL	5-061-C5-0	■			■			■			■		
0.5-10/20 µl	5-051-C5-0	■			■			■			■		
1-20/30 µl	5-116-C5-0	■	■		■	■		■	■		■	■	
1-50 µl	5-111-C5-0	■	■		■	■		■	■		■	■	
1-100 µl	5-112-C5-0	■	■		■	■		■	■		■	■	
1-200 µl	5-113-C5-0	■	■		■	■		■	■		■	■	

Filter Tips Compatibility
AHN Multi Channel Pipettes





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de