

COMPATIBILITY OF

AHN myTip Standard Tips

AHN myTip Macro Tips

AHN myTip Compatibility

with AHN Single Channel Pipettes

Tips Compatibility
 AHN Single Channel Pipettes

Pipettor		AHN																													
Brand		pipet4u [®] pro											pipet4u [®] Revolution																		
Model		variable/single channel					fix/single channel						variable/electronic single channel																		
Volume Range / Cat. No.		0.1-2.5 µl / 8-100-00-9	0.5-10 µl / 8-101-00-9	2-20 µl / 8-102-00-9	5-50 µl / 8-103-00-9	10-100 µl / 8-104-00-9	20-200 µl / 8-105-00-9	50-200 µl / 8-105-01-9	100-1000 µl / 8-106-00-9	200-1000 µl / 8-106-01-9	500-5000 µl / 8-107-00-9	2000-10000 µl / 8-108-00-9	5 µl / 8-001-00-9	10 µl / 8-002-00-9	20 µl / 8-003-00-9	25 µl / 8-004-00-9	50 µl / 8-007-00-9	100 µl / 8-008-00-9	200 µl / 8-009-00-9	250 µl / 8-010-00-9	500 µl / 8-011-00-9	1000 µl / 8-012-00-9	2000 µl / 8-013-00-9	2500 µl / 8-014-00-9	5000 µl / 8-015-00-9	10000 µl / 8-016-00-9	0.2-10 µl / 8-500-00-0	5-100 µl / 8-501-00-0	10-200 µl / 8-502-00-0	50-1000 µl / 8-503-00-0	100-5000 µl / 8-505-00-0
AHN myTip Standard Tips																															
Volume	Cat.No.																														
0.1-10 µl	1-000-81-0 1-000-C5-0 1-001-C5-0	■	■	■	■																					■	■	■			
0.1-10 µl XL	1-060-81-0 1-060-C5-0 1-061-C5-0	■	■	■	■																					■	■	■			
0.5-10/20 µl	1-050-81-0 1-050-C5-0 1-051-C5-0		■	■																											
1-200 µl	1-110-60-0 1-110-60-1 1-110-81-0 1-110-81-1 1-110-C5-0 1-110-C5-1 1-111-C5-0				■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
1-200 µl Gelfoader	1-130-70-0 1-130-C5-0 1-131-C5-0				■	■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■		
5-300 µl	1-103-81-0 1-103-C5-0 1-104-C5-0				■	■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■		
100-1000 µl short	1-209-80-2 1-209-C5-2 1-207-C5-2							■	■									■	■	■	■	■	■	■	■			■	■		
100-1000 µl long	1-200-80-0 1-200-80-2 1-200-C5-0 1-200-C5-2 1-201-C5-0							■	■									■	■	■	■	■	■	■	■	■	■	■	■		
AHN myTip Macro Tips																															
1-5 ml	1-210-84-0 1-214-84-0																														
1-5 ml	1-211-84-0 1-215-84-0																												■		
1-5.5 ml	1-212-84-0 1-216-84-0																														
1-10 ml	1-213-50-0 1-217-50-0											■														■					



AHN myTip Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN			
Brand		pipet4u® pro		pipet4u® Revolution	
Model	Volume Range / Cat. No.	variable/multi channel		variable/multi channel	
		8 Channel	12 Channel	8 Channel	12 Channel
		0.5-10 µl / 8-20100-9 5-50 µl / 8-202-00-9 30-300 µl / 8-203-00-9	0.5-10 µl / 8-401-00-9 5-50 µl / 8-402-00-9 30-300 µl / 8-403-00-9	0.2-10 µl / 8-520-00-0 2.5-50 µl / 8-521-00-0 15-300 µl / 8-522-00-0	0.2-10 µl / 8-530-00-0 2.5-50 µl / 8-531-00-0 15-300 µl / 8-532-00-0
AHN myTip Standard Tips					
Volume	Cat.No.				
0.1-10 µl	1-000-81-0	■	■	■	■
	1-000-C5-0	■	■	■	■
	1-001-C5-0	■	■	■	■
0.1-10 µl XL	1-060-81-0	■	■	■	■
	1-060-C5-0	■	■	■	■
	1-061-C5-0	■	■	■	■
0.5-10/20 µl	1-050-81-0	■	■	■	■
	1-050-C5-0	■	■	■	■
	1-051-C5-0	■	■	■	■
1-200 µl	1-110-60-0	■	■	■	■
	1-110-60-1	■	■	■	■
	1-110-81-0	■	■	■	■
	1-110-81-1	■	■	■	■
	1-110-C5-0	■	■	■	■
	1-110-C5-1	■	■	■	■
	1-111-C5-0	■	■	■	■
1-200 µl Gelloader	1-130-70-0	■	■	■	■
	1-130-C5-0	■	■	■	■
	1-131-C5-0	■	■	■	■
5-300 µl	1-103-81-0	■	■	■	■
	1-103-C5-0	■	■	■	■
	1-104-C5-0	■	■	■	■

Tips Compatibility
AHN Multi Channel Pipettes





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de